



Double Lumen Oocyte Collection Set

Designed for transvaginal aspiration of follicles to obtain oocytes. The set includes needle with tubing and bung for connection to a standard test tube. This needle allows continuous aspiration and flushing of the follicle thereby reducing operating time. All needles have echogenic markings for easy recognition under ultrasound.

This needle is available in 15-17 gauge, standard needle length 30cm. The standard length tubing set is 60cm or 85cm.

Materials: The needle and its connections are all manufactured from high grade surgical stainless steel. The tubing is manufactured from translucent PTFE. Both the needle and tubing are available in different lengths and gauges.

Sterilisation: Steam sterilised and designed for single use only.

Testing: All materials are tested to mouse embryo blastocyst survival to ensure a non-toxic, sterile product.

All products are manufactured in accordance with quality standards ISO 9001:2008, ISO 13485 and CE marking.