



## Embryo Transfer Catheter (Standard)

Our new range of catheters have been designed and manufactured in the United Kingdom to high industry standards of quality control.

It is reported that using soft embryo transfer catheters results in a significantly higher pregnancy rate as compared to firm catheters.

The range comprises of a soft catheter encased in an outer sheath. The catheter tip is rounded and smooth for non traumatic access. The outer sheath and catheter are marked with 1cm graduations to indicate degree of advancement.

The test catheter is available for use prior to embryo transfer and as a stylet only, for more complex procedures.

**Materials:** Catheter made from ultra soft non-embryo toxic polyurethane.

**Testing:** All materials are tested to mouse embryo blastocyst survival to ensure a non-toxic, sterile product.

**Sterilised:** Sterilised by using gamma irradiation and designed for single use only.

All products are manufactured in accordance with quality standards ISO 9001:2008, ISO 13485 and CE marking.

### Order details:

|                   |                      |                              |                            |
|-------------------|----------------------|------------------------------|----------------------------|
| Catalogue No:     | CAS/ET18<br>Catheter | CAS/ET18/TC<br>Test Catheter | CAS/ET18/ST<br>Stylet only |
| Quantity per box: | 10                   | 10                           | 10                         |
| Catheter length:  | 18cm                 | 18cm                         |                            |

Tel No.: +44 (0) 1372 720230 Fax No.: +44 (0) 1372 720260  
Email: sales@casmed.co.uk Website: www.casmed.co.uk