



**Toshiba  
Ultra**

**Siemens  
Short**

**Aloka  
Spring**

## Needle Guides

We hold a large range of metal reusable needles guides for use with ultrasound probes. They are designed to direct needles for accurate placement for oocyte collection.

### Features:

- Our instruments are designed and manufactured to a high specification in the UK and made of high grade surgical stainless steel
- Disposable sterile needles guides are also available for a wide range of ultrasound probes
- Guides accept needles from 18 gauge to 15 gauge
- Guides available for most ultrasound probes including Toshiba, Siemens, Aloka, etc
- Reusable guides reduce procedure costs
- Cleaning brushes available sharpness

**Materials:** High grade surgical stainless steel

**Sterilisation:** Produced non-sterile and re-usable. Can be sterilised by steam.

All products are manufactured in accordance with quality standards ISO 9001:2008, ISO 13485 and CE marking.

Tel No.: +44 (0) 1372 720230 Fax No.: +44 (0) 1372 720260  
Email: sales@casmed.co.uk Website: www.casmed.co.uk



**Aloka Spring  
mechanism**



**CASMED**  
INTERNATIONAL