



## Single Lumen Needle (Cyst Aspiration Needle)

Designed for transvaginal aspiration of cysts. This needle is available with or without a stylet. All needles have echogenic markings for easy recognition under ultrasound.

**Materials:** The needle and its connections are all manufactured from high grade surgical stainless steel.

**Sterilisation:** Steam sterilised and designed for single use only.

**Testing:** All materials are tested to mouse embryo blastocyst survival to ensure a non-toxic, sterile product.

Standard needle is 17 gauge x 31 cm. Also available in other gauges or lengths, please state when ordering.

All products are manufactured in accordance with quality standards ISO 9001:2008, ISO 13485 and CE marking.