



Veterinarian/Zoological Products

This double channel non-sterile needle may be used either laparoscopically or by ultrasound directed follicle aspiration for oocyte collection.

Features of this needle are:

Double lumen

- Needle designed and manufactured in the UK to high specification with proven record of quality control
- Marked near tip to enhance visualisation during ultrasound directed follicle aspiration
- Economical as the outer sheath can be easily replaced ensuring that the needle always has a sharp point without the expense of replacing the whole needle
- Tubing sets available separately at length required by customer. Easy to use as flushing and aspiration can occur without making adjustments to the device such as changing tubing
- Continuous flushing and aspiration reduces operating time
- Absence of "dead space" allows aspiration of complete follicle contents
- Needles and outers can be made from 10cm in length to 70cm – the length can be made to customer specification.

Concentric tube

- Combines optimum internal diameter for oocyte aspiration with acceptable outer diameter
- Needle is circular in cross section and there is no leakage if needle is rotated whilst in follicle

Needle point

- Specially ground in order to obtain maximum sharpness